

# **TUTORIALS to illustrate use of aXis2000 for the analysis of STXM \*.sdf data**

File: aXis2000-tutorials.doc

Author: Adam Hitchcock

Last changed: 11-Jan 2023

## **CONTENTS**

### **A. Image manipulation**

- A.1. Reading an ALS 532stxm image
- A.2. Converting transmission to optical density (OD)
- A.3. Determining profiles across an image
- A.4. Removing noise from images
  - A.4.1 Replacing lines
  - A.4.2 Smoothing images
  - A.4.3 Fourier filtering images
  - A.4.4 Clipping images
- A.5 displaying multiple Images

### **B Spectral Manipulation**

- B.1 Reading in a spectrum
  - B.1.1 Reading stxm5.3.2 multiple region spectra

### **C. Linescan Spectra processing**

- C.1 READ, Display, convert to OD
- C.2 Extracting and comparing spectra
- C.3 Fitting linescan spectra

- 7. Linescans~linefit 355-linescan-od.axb, the 3 components, select matrix, then pipa, then SAN spectra written out in step 6. Store the parameter file - 355linescan\_matrix\_pipa\_SAN.par
- 8. Spectra are plotted in buffers 1,2 , 3, which the distribution along the line for matrix, pipa and SAN are displayed in buffers 4, 5, 6. The residual of the fit is stored to disk (355-linescan-OD\_3SVD\_res.axb).

### **D. STACK processing**

#### **D. Stack processing**

[D.1 Read als5.3.2.2 stack & convert to OD](#)

[D.2 Extract reference spectra from the \\*od stack](#)

[D.3 Creating chemical maps from SVD fits](#)

-----

## A. Image manipulation

### A.1. READING AN ALS 532STXM IMAGE

```
* start aXis2000
* Read~STXM (sdf)
    Browse    - select ~\sdf-image.hdr
    OK
```

*Notes: The image is a 30x30 micron image of a wet cross-linked gel recorded at 396 eV. The upper right portion of the image is off-the-edge of the gel, passing through 2 silicon nitride windows and about 1 micron of water. (20320063)*

### A.2. CONVERTING TRANSMISSION TO OPTICAL DENSITY (OD)

[read in **sdf-image.hdf**] - tutorial #A.1.

- \* determine average signal in 'white' area (no sample) by either
  - \* move mouse over 'white area', read value from the cursor

-or-

\* **Images~Average pixels~region-** all pixels

*You should find a value of about 4500.*

- \* **Images~Convert\_to\_OD**

\* type 4500

The OD image is computed ( $OD = \ln(4500/\text{image})$ ); it is placed in buffer #0 – the default temporary buffer.

\* **Copy Buffer**, click on any empty thumbnail to copy the image to a permanent buffer.

### A.3. DETERMINING PROFILES ACROSS AN IMAGE

[read in **sdf-image.hdf**] - tutorial #A.1.

[convert to OD] – tutorial #A.2.

- \* **Images~profiles~linear**

- click at start of line, move mouse to end of line and click again
- select a buffer to store profile in; provide name if desired
- the OD intensity along the line is displayed

### A.4. REMOVING NOISE FROM IMAGES

#### A.4.1 Replacing lines

[read in **sdf-image-noisy.hdf**] - tutorial #1

- \* **Images~Replace line~manual-H**

- move cursor to just above the lower black horizontal line.

*The x, and y line number (lower left panel) should be 11*

- click mouse.

Line 11 is replaced by average of line 10 and line 12, and the result is moved to buffer #0.

- \* **Images~Replace line~manual-H**

- click on the upper part of the dark upper line – you will probably identify line 138.

- select 'no' then take the default in the next interaction box.

Which will replace line 138 with line 139.

- repeat to replace line 137 with average of lines 136 and 139.

*Note in this case, using the default (average adjacent lines) would leave a distortion in the image since two adjacent corrupt lines make up the upper dark line*

#### A.4.2 Smoothing images

[read in **image-to\_smooth.axb** using *Read~images~axis*]

- \* **Filter~smooth**

3- or 5-point smoothing reveals details of this image better.

#### A.4.3 Fourier filtering images

[read in **image-to\_clean.axb** using *Read~images~axis*]

- \* **Filter~clean (FT image filter)**

use cursor in rubber band mode (click at start, release, drag to end, click and release) to draw a rectangle top to bottom around the vertical line to the left of the central spot.

NB 2D Fourier transforms each quadrant is a replica of the other 3, so filtering out a pattern in one quadrant is all that is needed

*This is an image of a toner particle. The moiré pattern arises from beating of the frequency of vibrations of beamline optical element(s) against the ~kHz acquisition frequency.*

#### A.4.4 Clipping images

[read in **image-to\_clip.axb** using *Read~images~axis*]

- \* **Images~clip signal~histogram**

- a histogram of the z-values of this signal will appear. Select lower and upper bounds of data to be retained by moving mouse to lower limit of physically meaningful data (-15), click, then to upper limit (280), click.

- the pixels with values below the lower and above the upper limit are set to the closest limit and the result is placed in buffer #0.

- \* **Copy buffer** to move the result into a permanent buffer

*This is the map of styrene-co-acrylonitrile (SAN) from a fit a C 1s image sequence of a filled polyurethane. The haze is removed because pixels with large negative values are removed. These negative values from a fit are likely because the reference spectra are slightly wrong. **clip signal~histogram** removes these artificial results revealing more clearly the spatial distribution of the SAN particles.*

#### A.5 DISPLAYING MULTIPLE IMAGES

- \* **Read~images~axis image-to\_clean.axb** into buffer # 1

- \* Fourier filter this image to remove moiré noise (see Tutorial #A.4.4)

- \* **Copy buffer** to move the result into buffer # 2

- \* **Add~buffer #1** with -1 weighting to buffer #2

- \* **Copy buffer** to move the difference (removed noise) into buffer # 3

- \* view buffer # 1 by clicking on thumbnail #1

- \* **Filter~smooth** (3-point) the data; move to buffer #4

- \* **Display~thumbnails~4~common scale** selecting [1,2,3,4)

This type of comparison can help decide optimum processing. Note that the lineouts are active and give quantitative signal levels of all 4 displayed signals.

## B Spectral Manipulation

### B.1 READING IN A SPECTRUM

#### B.1.1 Reading stxm5.3.2 multiple region spectra

- start in buffer #1

- \* **read~STXM (sdf)** browse to select **11103106.hdr**

- in Region pull down menu, scroll to bottom, select all regions; OK  
Seven spectra will be read. Convert spectra in buffers 1 -6 to OD by:

- select buffer #1 **Spectra~Convert to~OD**; buffer 7 (Io)

- transfer result to buffer#1

- \* - select buffer #2 , 3 , 4, #5 (mtx), buffer #3 (pipa), #2 (san) and convert to OD using buffer #7 as Io. After generating each spectrum, transfer it to a suitable buffer to store for later spectral comparison

*This is a set of 7 point spectra acquired at same time on a 355 filler polyurethane sample. region 1: pipa (missed); 2, 6 - SAN, 3 - pipa; 4 - edge of SAN, 5 - matrix, 7, hole. (Stacks measured on this type of sample are used in the Advanced tutorial).*

The points at which these spectra were obtained can be viewed by:

- \* **read~images~graphics~png~image 355-point-spectra.png**

images recorded before and after these point spectra:

11103104 (285 eV) - highlights SAN, PIPA

11103105 (287 eV) - only SAN

11103107 (289 eV) - convert to OD; use Z-scale to check for damage

11103108 (285 eV) - should be same as 104

## C. Linescan Spectra processing

### C.1 READ, DISPLAY, CONVERT TO OD


1. Read~STXM (sdf) stxm532-linescan
2. extract Io (in region without sample)  
Linescans~Add lines~horizontal (y-values 1.0 – 1.5);  
Store in a permanent buffer (e.g. 6)
3. convert to OD, Linescans~normalize to Io buffer 6

### C.2 EXTRACTING AND COMPARING SPECTRA

4. Read~Images~axis 355-linescan-od.axb  
(This has already been converted to OD using the lower part of the line, off the polymer sample; it is the same sample as in D)
5. use linescans~Add lines~horizontal to extract the  
matrix (8.7 – 10.0), pipa (11.8-12.0) and SAN (13.3-14.00  
spectra (numbers are position along the line)
6. write each of these spectra to disk (Write~axis, and put name of  
the chemical species at start of header)

### C.3 Fitting linescan spectra

7. Linescans~linefit 355-linescan-od.axb, the 3 components, select  
matrix, then pipa, then SAN spectra written out in step 6. Store the  
parameter file - 355linescan\_matrix\_pipa\_SAN.par
8. Spectra are plotted in buffers 1,2 , 3, which the distribution along  
the line for matrix, pipa and SAN are displayed in buffers 4, 5, 6. The  
residual of the fit is stored to disk (355-linescan-OD\_3SVD\_res.axb).



## D. STACK processing

### D.1 Read als5.3.2.2 stack & convert to OD

532\_010817014 is a 38 energy C 1s stack measured on a polyurethane polymer consisting of 2 types of aromatic particles (SAN and pipa) in an aliphatic matrix [Ultramicroscopy 88 (2001) 33]

1. Convert STXM-sdf stack to an aXis2000 file  
- [Read~STXM \(sdf\)](#) \532\_010817014\10817014.hdr  
region#1 polymer region#2 - Io off microtomed section  
- output is \*.ncb, \*.dat', where \* is user selected name.  
the saved stack is then opened in [Stacks~Analyze~stack\\_process](#)
2. Align if needed (but this stack did not need alignment)
3. After reading region 2, select i: All and save the Io file
4. [Stacks~Analyze~stack\\_process](#) \*.ncb. Read Io from file; convert to OD; write out the stack \*.od.(ncb, \*,.dat)

### D.2 Extract reference spectra from the \*od stack

5. [Stacks~Analyze~stack\\_process](#) \*.od.ncb
6. Play back stack as a movie, Select I-region of largest particle. You should see a spectrum with 2 strong peaks, one at 285.1 eV ( $\pi^*(C=C)$ ) the other at 286.6 eV ( $\pi^*(C=N)$ ) save as \*-SAN
7. click on the spectrum display a the 286.6 eV peak. Some of the particles should disappear – they are the PIPA particles. Go back to 285.1 eV and select I: region of 1 or more of the small pipa particles. Save as \*-PIPA.
8. select a region with no particles – save as \*-matrix.

### D.3 Creating chemical maps from SVD fits

9. [Stacks~maps~SVD](#) - select \*.od.ncb
10. define 3 components – select the \*-matrix, \*-PIPA and \*-SAN as the reference spectra; save as matrix-pipa-san-14od.par
11. the Singular value decomposition fit will then execute, displaying the 3 reference spectra in buffers 1-3, and the 3 component maps in buffers 4-6, and the residual (data – fit) averaged over all photon energies in buffer 9.
12. [Images~clip~histogram](#) each of the component maps to remove -ve tails to improve contrast
12. [Display~RGB composite](#) (buffers 4, 5, 6) to make color composite of the 3 components. Select 'Yes' to rescale within each color

